

<b>Bauer Core Standard Protocol</b>		
Title: Hydroshear Protocol		
Pages: 2	Revision: 1.2	Date: 6/23/09
Author(s): Claire Reardon	Reviewers:	
Contact: <a href="mailto:claire@cgr.harvard.edu">claire@cgr.harvard.edu</a>	Comment:	

## 1. Purpose

This protocol describes how to use the GeneMachines (Genomic Solutions) Hydroshear to reproducibly fragment DNA. The protocol is designed as a reference and is not a substitute for training. Users must complete a training session before using any of the Bauer Core's instrumentation.

## 2. Materials

- 2.1. DNA sample(s) (see preparation instructions in section 4.1).
- 2.2. Wash Buffers (see preparation instructions in section 4.2).
- 2.3. Tubes to hold sheared DNA.
- 2.4. Shearing assembly.
  - 2.4.1. The Bauer Center provides standard shearing assemblies which can be used to shear DNA to sizes between 1 and 9KB. Users who would like to use either small (650B to 4KB) or large (4KB to 40KB) assemblies must provide their own.

## 3. Instrumentation

GeneMachines Hydroshear (now made by GenomicSolutions)

## 4. Reagent preparation

- 4.1. DNA sample(s)
  - 4.1.1. Resuspend the DNA sample in your buffer of choice (ddH<sub>2</sub>O or TE, pH 8.0).
    - 4.1.1.1. The concentration should be from 0.02ug/ul to 0.25 ug/ul.
    - 4.1.1.2. The volume should be from 40ul to 500ul.
 

Note that for sample volumes >300ul, a larger output tube must be used.
    - 4.1.1.2. The DNA **must** be completely in solution or it will clog the assembly.
 

Follow the instructions below to resuspend the DNA.
  - 4.1.2. Incubate for 30 minutes at 37°C, vortexing every 10 minutes.
  - 4.1.3. Immediately before shearing, centrifuge at 14,000 rpm for 15-25 minutes.
  - 4.1.4. If a pellet forms after centrifugation repeat steps 4.1.2 and 4.1.3.
  - 4.1.5. If a pellet still forms after repeated incubations, add more resuspension buffer.
  - 4.1.6. Whether or not you see a pellet, transfer the DNA solution to a clean tube, leaving behind the pellet if there is one.

## 4.2. Wash Buffers

- 4.2.1. Buffer 1: 0.2M HCl, 0.2 micron filtered.
- 4.2.2. Buffer 2: 0.2M NaOH, 0.2 micron filtered.
- 4.2.3. Buffer 3: ddH<sub>2</sub>O or TE, pH 8.0, 0.2 micron filtered.

## 5. Procedure

### 5.1 Setup.

- 5.1.1. Open the Hydroshear software.
- 5.1.2. Set shearing parameters.
  - 5.1.2.1. Volume: sample volume (40ul to 500ul).
  - 5.1.2.2. Number of Cycles: 20 cycles is standard.
  - 5.1.2.3. Speed Code: The standard assembly can shear at speeds from 3 to 40.  
Small speed codes indicate faster speeds which produce smaller fragments.
- 5.1.3. Click “Edit Wash Scheme” and choose how many times to wash in each buffer.  
A typical wash consists of 4 cycles in each buffer.

### 5.2. Shearing.

- 5.2.1. Hit “Start”.
- 5.2.2. Click “OK” to perform a wash unless the device has just been cleaned.
  - 5.2.2.1. Follow the prompts on the screen.
- 5.2.3. Follow the sample loading prompts until “turn the valve to output and click OK”.
  - 5.2.3.1. **Only turn the valve halfway so the handle is vertical**
  - 5.2.3.2. Click “OK” and watch the bubble rise as air is pushed out of the syringe.
  - 5.2.3.3. When the bubble rises to the top of the sample, turn the valve to output.
- 5.2.4. Continue following the prompts for shearing, sample ejection, and washing.
  - 5.2.4.1. Always perform a wash after each sample.

### 5.3. Shut down.

- 5.3.1. Always finish with a water rinse
- 5.3.2. Shut down the software.
- 5.3.3. Clean the area and remove any buffers or tubes that were used.

### 5.4 Tips and Tricks

- 5.4.1. To maximize sample recovery after shearing, aspirate and expel air instead of immediately performing a wash. The last bit of sample will be expelled.
- 5.4.2. When using the “standard” assembly, use a speed code of 3 or higher. When using a “small” assembly, determine the fastest possible speed (lowest number) by starting at 15 and trying faster and faster speeds. Using a speed that is too fast can cause leaks or over pressure failures.
- 5.4.3. When using a “small” assembly, use the manual settings to perform washes to allow for slower speed codes.
- 5.4.4. There are two ways to unclog the assembly:
  - 5.4.4.1 Perform the sonicating washes described in the manual.
  - 5.4.4.2. Turn the assembly around and run “backwards”.